



**Kinomatic Knee  
CT Protocol**

## TABLE OF CONTENTS

1	Introduction.....	3
2	Patient position .....	3
3	Image acquisition .....	4
	3.1 <i>Field of view (FOV)</i> .....	4
	3.2 <i>Image parameters</i> .....	5
	3.3 <i>Machine settings</i> .....	5
	3.4 <i>Scan ranges</i> .....	6
	3.5 <i>Image exportation</i> .....	7
4	Image delivery .....	7

## 1 INTRODUCTION

---

This document describes the guidelines that radiological centers should follow when performing CT-Tech Only scans required by Kinomatic LLC to use the Kinomatic Knee software.

The directives contained in this document are intended to improve the quality of the resulting CT scans. Radiology technicians are, in general, required to follow the instructions outlined in this document; variations to the settings and potential modifications to the protocol are to be discussed with and approved by Kinomatic LLC in advance.

Kinomatic LLC reserves the right to refuse scans performed with settings different from the suggested ones.

For more information, please do not hesitate to contact Kinomatic LLC at e-mail address: [radiology@kinomatic.com](mailto:radiology@kinomatic.com).

## 2 PATIENT POSITION

---

The patient must be in supine position at isocenter in the gantry. The legs must be in complete extension. No sponge or pillow should be placed beneath the knee or ankle.

The critical aspect is the position of the foot, which must be perpendicular to the table with the toe pointing straight up (see image below).



YES



NO

Importantly, the position of the foot should be secured to prevent motion and the consequent loss of accuracy to assess and correct any misalignment of the knee.

**TIP:** Tape the patient's feet to the table in the proper position to ensure they will not move during the scanning process and ensure the patient is comfortable and relaxed. This is critical for achieving a motionless scan.

### 3 IMAGE ACQUISITION

---

The acquisition consists of three (3) separate short **spiral axial scans**:  
**Scans must be the SAME FOV SIZE.**

- 1) Both hips
- 2) Both knees
- 3) Both ankles

Please note that scans of the hips and ankles are **required** to ensure an appropriate alignment of the legs. All three scans must be in the same coordinate system (frame of reference).

#### 3.1 FIELD OF VIEW [FOV]

---

Each acquisition must be centered and zoomed accurately to ensure the FOV maximizes the region of interest. Attention must be paid so that the outer bounds of the regions of interest are included in the FOV.

The whole of the **required bone regions** must be captured:

- **Hips:** femur heads to top of knees scan region (Left + Right)
- **Knees:** femurs, fibula and tibia contours (Left + Right)
- **Ankles:** medial and lateral malleoluses to bottom of knees scan region (Left + Right)

### 3.2 IMAGE PARAMETERS

Scans should be acquired in slices of minimum 512x512 pixels. The thickness of a single slice should be no more than 4mm for the hips and ankles and 1mm for the knees. The spacing between slices should be no larger than the slice thickness; a slight overlapping is allowed. UTILIZE BONE WINDOWS TO MAXIMIZE VIEW OF BONES

The following table summarizes the recommended image settings.

**NOTE: ALL 3 IMAGES MUST BE OF THE SAME SIZE**

	KNEES	ANKLES & HIPS
ROWS & COLUMNS	512 x 512 pixels	512 x 512 pixels
SLICE THICKNESS	0.5 - 1mm	1 - 4mm
SPACING BETWEEN SLICES	0.5 - 1mm	1 - 4mm
FOV	Min 350mm/Max 500mm	Min 350mm/Max 500mm

The FOV should be as small as possible, as long as the articulation is completely demonstrated.

### 3.3 MACHINE SETTINGS

The following settings are intended to maximize the quality of the images. If different values are required on your machine, please inform Kinomatic LLC before proceeding with the acquisition.

<b>KiloVolt Peak (KVP)</b>	120 KV or higher
<b>X-Ray Tube Current</b>	120 mA or higher

### 3.4 SCAN RANGES

---

The following figure depicts the scans required of both sides:



#### HIPS

The whole femoral head and continue to **15 cm** below the lesser trochanter

#### KNEES

Both femoral and tibial parts of the knee joint must be included in the knee acquisition. The scan must be extended to at least **15 cm** towards the hip and **10 cm** towards the ankle, past the bottom of the tibial tuberosity.

#### ANKLES

The ankles scan must include at least **15 cm** of the tibiae and extended past the calcaneus

### 3.5 IMAGE EXPORTATION

---

The three required scans of the hips, knees and ankles must be **exported as original scans**. In particular, derived (re-sampled) series will be rejected. The DICOM dataset provided should contain only the series of interest to minimize confusion and likelihood of error. Regarding the Reconstruction Convolution Kernel, we suggest using “bone kernel” for all three scans (hips, knees and ankles).

Both the study and the series should be named according to what they represent. As a rule of thumb, include the patient’s name and affected side in the study description, and the affected side and part in the description of each series, e.g.:

Study: “John Smith, knee study”

Series 1: “Hip”

Series 2: “Knee”

Series 3: “Ankle”

In general, any localizers automatically produced by the machine should be excluded from the media. If this is not possible, please clearly distinguish between localizers and actual scans by providing a series consisting of localizers only. **DO NOT include localizers within the series.**

The image data exported should be in uncompressed RAW format (MONOCHROME2) or compressed using lossless algorithms (i.e. JPEG2000 lossless).

## 4 IMAGE DELIVERY

---

The DICOM dataset of a Kinomatic CT case can be uploaded to an online portal by the radiological center.

Alternatively, the DICOM dataset can be stored on a CD or DVD. A label must be applied including the patient’s name, affected side, acquisition date and, if applicable, planned surgery date.

The storage media can then be mailed to:

Kinomatic LLC  
136 W. Branch St.  
Arroyo Grande, CA 93420

Do not hesitate to contact Kinomatic LLC for further assistance at: [radiology@kinomatic.com](mailto:radiology@kinomatic.com).





REVOLUTIONARY TECHNOLOGY  
SUPERIOR OUTCOMES

[KINOMATIC.COM](http://KINOMATIC.COM)



**Kinomatic LLC**  
136 West Branch Street  
Arroyo Grande, California 93420  
Phone (800) 420-1701  
[info@kinomatic.com](mailto:info@kinomatic.com)

All trademarks and registered trademarks are the property of their respective owners.  
This document is not intended for use outside of the US market.  
Please verify approval of the devices described in this document with a Kinomatic LLC representative.

Rev. 2024-01