

	Knee VIP™ and Hip VIP™ Instructions for Use
---	--

INTENDED USE

Knee VIP™ and Hip VIP™ are preoperative surgical planning software. They are intended to assist qualified orthopedic surgeons in the preoperative planning of orthopedic surgical procedures.

INDICATIONS FOR USE

Knee VIP™ and Hip VIP™ are pre-operative planning tools for orthopedic surgery. The plans are created from patient pre-dimensioned CT scans from PACS-systems, utilizing implant specific surgical techniques with surgeon preferences. The Knee VIP™ and Hip VIP™ system contain an integrated database of orthopedic hip and knee implant geometries that are overlaid to assist surgeons in their planning of orthopedic hip and knee surgeries. The software application enables the healthcare professional to view the patient’s customized preoperative plan by means of a Virtual Reality (VR) headset. Finalized plans can be viewed with VR or printed to a PDF report. Clinical judgment and experience with the software are required for its successful use.

CONTRAINDICATIONS

Knee VIP™ and Hip VIP™ plans are not intended for patients who are not physically able to receive a CT scan.

USE

Detailed usage is outlined in the Knee VIP™ and Hip VIP™ User Manual.











REUSE


Knee VIP™ and Hip VIP™ are custom surgical plans meant for one time use only. They are made for one patient and one procedure.

CAUTION

Federal Law (USA) restricts this device to sale by or on the order of a physician.

PRODUCT LABELS LEGEND

				
Catalog Number	Lot Number	Serial Number	UDI Carrier	Caution required when using device: Refer to IFU
	 eIFU Indicator			
For Prescription Use Only	Consult eIFU using given URL before using device	Expiration Date	Manufacturer	Manufacturing Date

Kinomatic LLC

 136 W Branch St
 Arroyo Grande, CA 93420 U.S.A.
 Phone: (800) 420-1701
 Website: <https://kinomatic.com/>

Paper IFU will be mailed upon request